

Bacterial milk testing for dairy plants who want:

- Insight into milk quality.
- Critical thermodurics detected.
- 3,5-hour sample-to-result.
- Improved payment schemes.
- Optimal usage of their milk.

Milk quality is one of the most crucial factors in dairy production. Insight into the critical bacteria in your milk means not overpaying suppliers and helps to determine which milk is best for specific products. Lacti Range provides qPCR tests to identify live bacteria in milk before and after heat treatment. We know that efficiency is essential, and with rapid test results in only 3,5 hours, our tests minimize any holdup in production and ultimately save money.

For more info visit www.dna-diagnostic.com or contact Tel +45 87 32 30 50 · info@dna-diagnostic.com

